



PARTNERS
FOR ENDOSCOPY

Triozyme Enzymatic Cleaning Solution



Current guidelines support the use of enzymatic detergent in cleaning by the patient bedside – during the manual clean and in an AER.

Our enzymatic detergent can be used in all stages of the endoscope cleaning process and can even be used in ultrasonic baths. Triozyme is a low foaming, **non-irritating** enzymatic effective over a broad temperature spectrum and demonstrated to provide superior efficacy in the removal of soil in tests against other enzymatic cleaners, without the need for elevated temperature.

(Available from Partners in 5 litre bottles complete with a pelican pump screw on top - an automatic dosing pump can be supplied free of charge in some instances)

Product Code: 40003

Box of 4 x 5litres

For more information please call **01782 594164**
or via email **sales@pfmedical.co.uk**

Hartley House, Galveston Grove, Fenton, Stoke-on-Trent ST4 3PE www.pfmedical.co.uk

All prices are valid for 30 days and exclude VAT and Carriage

Last updated June 2018